

# **ACTIVE VIBRATION ISOLATION SYSTEM INSTALLATION REPORT**

- 1. Measurement Details**
- 2. Equipment Information**
- 3. Vibration Isolation System**
- 4. Installation Photo**
- 5. Summary**
- 6. Results**

# 1. Measurement Details

---

- **Measurement Date**

January 13, 2020

- **Measurement Devices**

1. LAN-XI Data Acquisition Hardware

- Brüel & Kjær 3050-B-040

2. Data Analysis Software

- Brüel & Kjær PULSE LAB SHOP 14

3. Sensors

- PCB Accelerometer

- Model: 393B05

- **Location**

1<sup>st</sup> Floor

- **Irregular Event**

The installation site is located near a subway station within 3 m.

## 2. Equipment Information

- **Manufacturer**

ZEISS

- **Model**

Gemin FESEM

- **Vibration Specification**

Allowable horizontal vibration values

up to 10 Hz : less than 5  $\mu\text{m/s}$  rms (VC-E)

10 – 60 Hz : less than 10  $\mu\text{m/s}$  rms (VC-D)

above 60 Hz : less than 14  $\mu\text{m/s}$  rms (VC-C)

Allowable vertical vibration values

up to 10 Hz : less than 4  $\mu\text{m/s}$  rms (VC-E)

10 – 60 Hz : less than 14  $\mu\text{m/s}$  rms (VC-C)

above 60 Hz : less than 20  $\mu\text{m/s}$  rms (VC-C)


### 3. Vibration Isolation System

**Model: DVIA-MB1000**



Platform Dimensions	1140 x 910 x 224 mm	
Load Capacity	500 - 1700 kg	
Actuator	Electromagnetic Actuator	
Maximum Actuator Force	Vertical: 40N, Horizontal: 20 N	
Degrees of Freedom	6 degrees	
Active Isolation Range	0.5 - 100 Hz	
Vibration Isolation at 2 Hz	≥90%	
Vibration Isolation at 10 Hz	≥99%	
Input Voltage (V)	AC100 - 240V / 50 - 60 Hz / 1A	
Power Consumption (W)	Maximum 110W, <50 W in normal operation	
Operating Range	Temperature (°C)	5 - 50 °C
	Humidity (%)	20 - 90%
Required Air Pressure	≥ 0.5 MPa (5 bar)	


## 4. Installation Photo



## 5. Summary

The installation site is located near a subway station, caused the low frequency floor vibration in vertical and horizontal axis that the measured floor vibration levels did not meet the ZEISS Gemini SEM vibration specification. DAEIL SYSTEMS's active vibration isolation system reduced the floor vibrations in Z, X and Y axis, succeeded in meeting the vibration specification.


## 6. Results – Z axis (Vertical)



The subway trains caused a significant amount of vibration in Z axis at 25 Hz, at which point the floor vibration exceeded the allowable vertical vibration values for ZEISS Gemini SEM. Our active vibration isolation system reduced the 25 Hz floor vibration, meeting the vertical vibration specification.




## 6. Results – X axis (Left to Right)



The subway trains caused floor vibration at 2 Hz and 25 Hz in x axis, which exceeded the horizontal allowable vibration values. Our active vibration isolation system reduced the floor vibration from VC-D to VC-G, meeting the electron microscope's vibration specification.

## 6. Results – Y axis (Front to Back)



The subway trains caused floor vibration at 2.25 Hz in y axis, which exceeded the horizontal allowable vibration values. Our active vibration isolation system reduced the floor vibration from VC-D to VC-G, meeting the electron microscope's vibration specification.